




## Introducer Set



- Dilator locked clockwise does not move backward and allows easy operation.
- Close contact between the tip of sheath and dilator developed through our long term of studies and experience allows safe and smooth insertion into vessel.
- Supercath Seldinger Needle is manufactured by one piece molding process developed by us ensures easy guidewire insertion.

**Packing :** Box of 10 pcs.

**Size :** Fr 6, 7, 8.

### LIFELINE SYSTEMS PVT. LTD.

Regd. Office : A.D. House, J-161/1, Gautam Nagar, New Delhi-110049. INDIA.  
Sales Office : F-63, Okhla Industrial Area Phase-1, New Delhi-110020. INDIA.  
Phone : +91-11-41640924 - 6, 40523028  
Fax : +91-11-40523027, 26810652  
Email : info@lifelinedelhi.com  
URL : www.lifelinedelhi.com

Branch Offices : 111/506, Navjivan Society, Lamington Road, Mumbai-400008. India.  
Phone : +91-22-23096796  
Sundaram Building, 1-D, 1st Floor, 46 Rafi Ahmed Kidwai Road,  
Kolkatta-700016. India. Phone : +91-33-22299932